UCI Beall **Applied Innovation**

Research Translation Group

Request Information

Research Translation Group

Available Technologies

Contact Us

Permalink

A Combined Microfluidic and Fluorescence Lifetime Imaging(FLIM) Platform to Identify Mammalian Circulating Cancer Cells in Whole Blood

Tech ID: 28924 / UC Case 2017-815-0

BRIEF DESCRIPTION

Separating and classifying circulating cancer cells from whole blood using a single cell trap microfluidic platform coupled with label free fluorescence life time imaging.

FULL DESCRIPTION

Blood-borne metastasis caused by circulating cancer cells (CTC) is the cause of 90% of cancer related deaths. It has been shown that circulating cancer cells can be detected and classified by a liquid blood biopsy. Liquid biopsy utilizes microfluidics and imaging as a way to detect single CTCs from dense whole blood. Unfortunately, these current microfluidic biopsy techniques lack the sensitivity to detect CTCs in whole blood due to non-selective cell sorting nature of these methods.

Researchers at the University of California have developed a robust sensitive lab-on-a-chip platform to detect CTCs through a combined microfluidics and label-free approach via fluorescence lifetime imaging microscopy (FLIM). FLIM is a non-invasive methodology that can differentiate innate biological fluorescent absorbers based on the time it takes for the absorber's fluorescence brightness to fall. A FLIM-Phasor plot that contains all fluorescent lifetime(s) of a sample can provide a straight forward visualization of the different biochemical compounds and their chemical state in the sample. Utilizing such a FLIM-Phasor plot has shown clear differentiation between different cancer cells providing greater diagnostic information regarding the isolated CTCs.

SUGGESTED USES

- Separating circulating tumor cells from whole blood
- » Diagnosing sickle cell anemia
- » Differentiate cancer and leukemia cells

ADVANTAGES

- » Revealing molecular-level signature of tumor
- » Low cost, disposable, ease of mass production
- » Multifunctional diagnostic platform for different blood related diseases

Туре

PATENT STATUS

Country

Number

Dated

Case

CONTACT

Alvin Viray aviray@uci.edu tel: 949-824-3104.



OTHER INFORMATION

CATEGORIZED AS

- >> Optics and Photonics
 - All Optics and **Photonics**
- » Medical
 - » Devices
 - » Diagnostics
 - Disease: Blood and Lymphatic System
 - >> Imaging
 - >>> Research Tools
 - >>> Screening

RELATED CASES

2017-815-0

STATE OF DEVELOPMENT

Invention has completed the conceptualization phase, simulation and preliminary experimental stage. Currently, the invention is undergoing design optimization to improve cell entrapment.



5270 California Avenue / Irvine,CA 92697-7700 / Tel: 949.824.2683



© 2017 - 2023, The Regents of the University of California Terms of use Privacy Notice