

Technology Development Group

Available Technologies

Contact Our Team

CONTACT

Permalink

Request Information

Guided Magnetic Nanospears For Targeted And High-Throughput Intracellular Delivery

Tech ID: 29468 / UC Case 2018-436-0

SUMMARY

UCLA researchers in the Department of Chemistry & Biochemistry and Department of Molecular & Medical Pharmacology have developed novel magnetic nanostructures that can be used to carry and/or deliver biomolecular cargo intracellularly to cells.

BACKGROUND

Current approaches for high-throughput and targeted intracellular delivery of biomolecules involve the use of viruses, external electric fields, or harsh chemical reagents, which are either costly, inefficient, or apply undesirable stresses or toxicities to the cells. Improving the safety, speed, cost effectiveness, and efficiency of intracellular delivery remains a challenge in the field of cell and molecular biology, and gene editing systems for clinical applications.

Membrane disruption-based approaches using sharp, needle-like nanostructures can physically penetrate flexible cell membranes to deliver biomolecules to cells efficiently, with minimal impact on cell viability and metabolism. Most nanoneedle platforms consist of arrays grown on planar substrates that support the growth of adherent cells, but they have problems releasing the modified cells after transfection and collecting them for further studies. Recently developed nano-/micromotor systems can be internally or externally powered to move in liquid environments, but these nanosystems have limited precision in terms of their guidance and biocompatibility due to byproducts from catalytic reactions that propel the nano-/micromotor structures to their targets.

INNOVATION

Researchers at UCLA have developed magnetic nanostructures that can be configured readily for single-cell modification, or they can be scaled progressively for direct and highly efficient high-throughput intracellular delivery. These biocompatible nanomaterials can be guided precisely to target cells without the need for chemical propellants *via* manipulation of locally applied magnetic fields. The guided nanospears can be manufactured and deployed in large scales while achieving exceptional transfection efficiencies, maintaining cell viability for applications in stem cell biology and the development of next-generation gene and cellular immunotherapies.

APPLICATIONS

> Carry and deliver biomolecular cargo (i.e. DNA plasmids, siRNA, proteins, targeted endonucleases) intracellularly to cells

ADVANTAGES

The position, orientation, and speed of the magnetic nanospears can be manipulated remotely via an external magnetic field, without needing chemical reagent

- > The surfaces of the nanospears can be easily configured to carry a variety of biomolecular cargos
- The nanospears can be used to selectively transform neighboring cells with single-cell precision, or to transform a large population of cells

STATE OF DEVELOPMENT

The described guided nanospears have been used to successfully carry DNA plasmids into target cells.

Group ncd@tdg.ucla.edu tel: 310.794.0558. INTRODUCING UC TechAlerts New technology matches delivered to your email at your preferred schedule

UCLA Technology Development

your email at your preferred schedule SEARCH > & SAVE SEARCH Learn More

INVENTORS

Weiss, Paul S.

OTHER INFORMATION

KEYWORDS Nanosubstrate-mediated delivery, nanomaterial, intracellular delivery, transformation, transfection

CATEGORIZED AS

Materials & Chemicals

- Nanomaterials
- Medical
 - Research Tools
- Nanotechnology
 - NanoBio

RELATED CASES 2018-436-0

Country	Туре	Number	Dated	Case
United States Of America	Published Application	20200347409	11/05/2020	2018-436

ADDITIONAL TECHNOLOGIES BY THESE INVENTORS

- Determining Oil Well Connectivity Using Nanoparticles
- Multiple-Patterning Nanosphere Lithography
- ▶ High-Throughput Intracellular Delivery of Biomolecular Cargos via Vibrational Cell Deformability within Microchannels
- Scalable Lipid Bilayer Microfluidics for High-Throughput Gene Editing
- ▶ High-Throughput Microfluidic Gene-Editing via Cell Deformability within Microchannels

Gateway to Innovation, Research and Entrepreneurship

UCLA Technology Development Group

© 2018 - 2020, The Regents of the University of California



https://tdg.ucla.edu Tel: 310.794.0558 | Fax: 310.794.0638 | ncd@tdg.ucla.edu

10889 Wilshire Blvd., Suite 920,Los Angeles,CA 90095

Terms of use Privacy Notice